

Angiostamp™800

NIR molecular imaging agent for oncology


Angiostamp™800 is a targeted contrast agent for fluorescence and photoacoustic imaging of tumors.

Angiostamp™800 is an organic NIR imaging agent targeting the $\alpha v \beta 3$ integrin, a cell surface receptor involved in cell adhesion and migration. The $\alpha v \beta 3$ integrin is overexpressed in particular on endothelial cells during neoangiogenesis as well as on many tumor cells types.


4 cyclic RGD for $\alpha v \beta 3$ recognition

Decapeptide platform


NIR absorbing DYE



Angiostamp™800



Photoacoustic absorption spectrum



Fluorescence absorption and emission spectra

Angiostamp™800 for Photoacoustic imaging in oncology


Orthotopic glioblastoma model / multimodal investigation using Angiostamp™800

Bioluminescence




→ Follow-up of cancer cells proliferation

2D PAI / US / μ CT



→ Monitoring of Tumor specific biomarker expression in a detailed anatomical context

3D PAI / US / μ CT



Angiostamp™800 was first developed for *in vivo* tumor detection by fluorescence imaging [1]. Since then, it has been widely characterized and was shown to provide high tumor contrast in various animal models including lung adenocarcinoma [2, 3], ovarian carcinoma [4, 5], peritoneal carcinomatosis [6, 7], fibrosarcoma [8], osteosarcoma [9], head and neck squamous cell carcinoma [10] and bone metastasis from breast cancer origin [11].


Angiostamp™800

NIR molecular imaging agent for oncology

Angiostamp™800 for Photoacoustic imaging in oncology


Orthotopic prostate model / multimodal investigation using Angiostamp™800

Bioluminescence




→ Follow-up of cancer cells proliferation

2D PAI / US / μ CT




→ Monitoring of Tumor specific biomarker expression in a detailed anatomical context

3D PAI / US / μ CT




Orthotopic bone metastases model / multimodal investigation using Angiostamp™800

Bioluminescence



→ Follow-up of cancer cells proliferation

2D PAI / US / μ CT




→ Monitoring of Tumor specific biomarker expression in a detailed anatomical context

3D PAI / US / μ CT




Angiostamp™800 can increase up to 43 times the tumor contrast depending on $\alpha\beta 3$ expression level.



1. Garanger, E., et al., *New multifunctional molecular conjugate vector for targeting, imaging, and therapy of tumors*. Mol Ther, 2005. **12**(6): p. 1168-75.
2. Koenig, A., et al., *In vivo mice lung tumor follow-up with fluorescence diffuse optical tomography*. J Biomed Opt, 2008. **13**(1): p. 011008.
3. Koenig, A., et al., *Fluorescence diffuse optical tomography for free-space and multicolor studies*. J Biomed Opt, 2010. **15**(1): p. 016016.
4. Jin, Z.H., et al., *Noninvasive optical imaging of ovarian metastases using Cy5-labeled RAFT-c(-RGDfK)-4*. Mol Imaging, 2006. **5**(3): p. 188-97.
5. Mery, E., et al., *Intraoperative fluorescence imaging of peritoneal dissemination of ovarian carcinomas. A preclinical study*. Gynecol Oncol, 2011. **122**(1): p. 155-62.
6. Keramidas, M., et al., *Intraoperative near-infrared image-guided surgery for peritoneal carcinomatosis in a preclinical experimental model*. Br J Surg, 2010. **97**(5): p. 737-43.
7. Josserand, V., et al., *Electrochemotherapy guided by intraoperative fluorescence imaging for the treatment of inoperable peritoneal micro-metastases*. J Control Release, 2016. **233**: p. 81-7.
8. Wenk, C.H., et al., *Near-infrared optical guided surgery of highly infiltrative fibrosarcomas in cats using an anti- α v β 3 integrin molecular probe*. Cancer Lett, 2013. **334**(2): p. 188-95.
9. Dutour, A., et al., *Targeted imaging of α (v) β (3) expressing sarcoma tumor cells in vivo in pre-operative setting using near infrared: a potential tool to reduce incomplete surgical resection*. Bone, 2014. **62**: p. 71-8.
10. Atallah, I., et al., *Near-infrared fluorescence imaging-guided surgery improves recurrence-free survival rate in novel orthotopic animal model of head and neck squamous cell carcinoma*. Head Neck, 2016. **38 Suppl 1**: p. E246-55.
11. Bellanger, A., et al., *The critical role of the ZNF217 oncogene in promoting breast cancer metastasis to the bone*. J Pathol, 2017. **242**(1): p. 73-89.

Contact information :

Optimal Grenoble
Rond point de la Chantourne
38700 La Tronche
FRANCE

 **+33 (0)4 76 54 94 10**

 **contact@imagerie-optimal.fr**