

CYBERSECURITY ALERT

Dear VisualSonics customer,

FUJIFILM VisualSonics, Inc. ("VSI") is aware of the ZombieLoad and ZombieLoad V2 vulnerabilities in Intel CPUs, and has evaluated the potential impact on its imaging systems. Intel, manufacturer of the processors used in VSI imaging systems, has communicated that certain processor models are vulnerable to ZombieLoad and ZombieLoad V2. After a thorough investigation, VSI can confirm that the following currently supported pre-clinical imaging systems contain a potentially vulnerable Intel processor:

Vevo 3100	Vevo LAZR-X	Vevo 2100
Vevo LAZR	Vevo 1100	

To date, there are no known exploits (other than laboratory experiments) that can take advantage of ZombieLoad and ZombieLoad V2. In order to exploit the vulnerabilities, an attacker must be able to run specially crafted code directly on the system. VSI imaging systems are closed devices which do not allow users to run other applications or code. Although the underlying CPU and operating system combination may be affected by these vulnerabilities, the fact that the systems are closed means that they are not vulnerable - there is no vector available through which to exploit the systems.

VSI is committed to ensuring that its products are safe, secure, and reliable. When stable security patches are released by Microsoft and / or Intel, VSI will thoroughly test to ensure safety and then integrate those patches for the following products:

Vevo 3100	Vevo LAZR-X	Vevo 2100
Vevo LAZR	Vevo 1100	

VSI has determined that the Vevo 770 imaging system cannot be patched.

VSI will notify customers when patched versions of imaging system software are released. Should you have any questions in the interim, please contact VSI technical support:

Web: www.visualsonics.com

U.S.A. or Canada:

+1 877.657.8118 or

+1 888.482.9449

Outside U.S.A. or Canada:

+1 425.951.1330

or call your local representative

Fax: +1 425.951.6700

Asia Service Center: +65 6380 5581

E-mail: support@visualsonics.com

European Service Center:

Main: +31 20 751 2020

English support: +44 14 6234 1151

French support: +33 1 8288 0702

German support: +49 69 8088 4030

Italian support: +39 02 9475 3655

Spanish support: +34 91 123 8451